

INSTRUCTIONS

Chrometra Tubulin ExM

Polymerizable Taxane derivatives

Product Numbers	Description
ER-15 range	Reagent for staining of tubulin structures in microscopy. Contents: Reactive taxane derivative (100 or 300 reactions) Storage Upon receipt store product at -20°C.

Introduction

Chrometra Tubulin ExM for staining of tubulin networks in polymeric gels and general fluorescence microscopy.

The Chrometra Tubulin ExM kit is used to introduce reactive labels onto microtubules. By incorporating a polymerizable moiety into the probe, the conjugate becomes covalently anchored to the swellable gel during polymerization, preserving the tubulin signal through digestion and isotropic expansion. The reagent is available containing either a reactive azide or a biotin.

Kit Contents

- 1 or 3 vials of Tubulin ExM grafting reagent (100 or 300 reactions, 30 reactions for samples).
- 1 or 3 vials of fluorescent Post Stain reagent
- Optional: 20 mM N-Ethyl Maleimide (NEM) (see recipe below).

Reagents not included but which may be required:

- Glutaraldehyde
- NaBH₄ as a 0.15% (w/v) solution in PBS
- PBS Buffer
- PBS Buffer containing 0.1% Tween-20
- PBS + 1% BSA
- PEM solution (80 mM PIPES pH 6.8, 5 mM Ethylene glycol tetraacetic acid (EGTA), 1 mM MgCl₂ and 0.3% Triton X-100 (see recipe below).

Reagent preparation

The grafting reagent is supplied dry. Reconstitute it in DMSO or methanol by adding 100 µL of solvent to the tube and incubating at 37 °C with gentle shaking for 5 min. Immediately before staining, dilute 1 µL of this stock into 50 µL of an appropriate buffer (e.g., PBS; not provided). Reconstitute the Post stain reagent in DMSO by adding 100 µL of solvent to the tube and incubating at 37 °C with gentle shaking for 5 min. 1 µL of this stock corresponds to one unit. Store both reagents frozen.

Procedure for tubulin grafting (default)

1. Incubate the sample in 3 µM docetaxel probe for 1 h at 37 °C.
2. Fix in 0.5% GA in PEM + 0.3% Triton X-100 for 10 min.
3. Rinse once with PBS.
4. Quench residual GA with 0.15% NaBH₄ in PBS for 7 min at RT.
5. Wash the cells with PBS buffer (3 times, for 5 minutes)
6. Optional: Incubate in 20 mM NEM (PBS) for 30 min to block free thiols
7. Wash the cells with PBS buffer (3 times, for 5 minutes)
8. Proceed to staining through Click Labeling (Procedure A, B) or direct streptavidin staining (Procedure C)

Labeling Procedure A: Click Labeling (Biotin/Streptavidin example, for use with Tubulin ExM Biotin)

1. Incubate with 10 μ M of in PBS + 1% BSA.
2. Wash with PBS + 0.1% Tween-20. (3 times, for 5 minutes)
3. Incubate with fluorescent streptavidin (15 μ g/mL) in PBS + 1% BSA for 1 h.
4. Wash the cells with PBS buffer (3 times, for 5 minutes)
5. Optional: Nuclear counterstain: Incubate with 1 μ g/mL DAPI for 1 min.
6. Wash the cells with PBS buffer (3 times, for 5 minutes)
7. Proceed to gelation or imaging

Labeling Procedure B: Click Labeling (Dye example, for use with Tubulin ExM Post)

1. After the default tubulin grafting, incubate with 1 unit of Post Stain reagent dye in PBS + 1% BSA.
2. Wash the cells with PBS buffer + 0.1% Tween-20 (3 times, for 5 minutes)
3. Proceed to gelation or imaging

Labeling procedure C: Direct streptavidin staining (for use with Tubulin ExM Biotin)

1. Incubate with fluorescent streptavidin (15 μ g/mL) in PBS + 1% BSA for 1 h.
2. Wash the cells with PBS buffer (3 times, for 5 minutes)
3. Proceed to gelation or imaging

Example Calculations for stock solutions:

0.15% (w/v) NaBH₄ in PBS:

Weigh 1.5 mg NaBH₄ and bring up to 1.00 mL with PBS (prepare fresh, cold if possible).

20 mM N-ethylmaleimide (NEM) in PBS:

Bring the provided tube of NEM up to 1.00 mL with PBS.

PEM solution (10 mL recipe):

1. Weigh the following reagents:
 - a. PIPES (80 mM; 0.8 mmol): PIPES free acid (MW 302.37): 242 mg, or PIPES disodium (MW 366.35): 293.1 mg.
 - b. EGTA (5 mM; 50 μ mol): EGTA free acid (MW 380.35): 19.0 mg or EGTA disodium (MW 436.31): 21.8 mg.
 - c. MgCl₂ (1 mM; 10 μ mol): Anhydrous (MW 95.21): 0.952 mg, or Hexahydrate (MW 203.30): 2.03 mg. (Practical alternative: add 10 μ L of a 1 M MgCl₂ stock.)
 - d. Triton X-100 (0.3% v/v): 30 μ L
2. Prepare the solution:
 - a. Dissolve PIPES and EGTA in ~8 mL ultrapure H₂O.
 - b. Adjust to pH 6.8 with NaOH.
 - c. Add MgCl₂ (weighed or from 1M stock).
 - d. Add 30 μ L Triton X-100, mix gently.
 - e. Bring volume to 10 mL with H₂O
 - f. Optional: filter-sterilize; store cold.

Example gelation procedure (4x Expansion)

1. Incubate samples in 0.1 mg/ml crosslinker (e.g.; Acrylink NHS (0.35 mM), Acrylink TFP (0.3 mM), AcX (0.35 mM)) in PBS for 90 min at RT.
2. Wash the cells with PBS buffer (3 times, for 5 minutes)
3. Prepare gelation mix (monomer + initiators): 2 M NaCl, 8.625% (w/w) sodium acrylate, 2.5% (w/w) acrylamide, 0.15% (w/w) N,N'-methylenebisacrylamide, 0.01% 4-hydroxy-TEMPO.
4. Immediately before use, add 0.15% (w/w) TEMED and 0.15% (w/w) APS to the required aliquot. Mix gently.
5. Quickly rinse each sample with 50 μ L of the freshly initiated monomer solution to pre-infuse.
6. Assemble the sample in a gelation chamber.

7. Cover the sample in gelation solution
8. Polymerize in a humidified incubator at 37 °C for 90 min (avoid evaporation).
9. Digest by incubating the gels with Proteinase K, 8 U/mL in digestion buffer (50 mM Tris, pH 8.0; 1 mM EDTA; 0.5% Triton X-100; 0.8 M guanidine HCl) at 37 °C for 3 h.
10. Optional nuclear counterstain: Stain with DAPI, 1 µg/mL in PBS for 1 min; Wash the cells with PBS buffer (3 times)
11. Transfer gels to deionized H₂O for 10 min.
12. Replace with fresh H₂O and repeat 4× (or until dimension change plateaus) to reach full 4× expansion.

Storage and Safety

The Tubulin ExM reagent can be stored at 4°C for 3 months. For longer storage the conjugate can be stored at -20°C. The conjugate should always be stored in the dark.

Tips and troubleshooting

- To avoid photobleaching of dilute dyes, minimize the exposure of fluorescently labeled specimens to light.
- Dye concentration: When following the standard guidelines, one unit of dye reagent in a staining buffer solution of 100 µl will result in a 5 µM concentration. Dye working concentrations can be reduced to accommodate for larger volumes of staining solution. In such case, extend incubation to 2 hours at 37°C.
- The addition of 0.1 mg/ml general linker (e.g.; Acrylink NHS (0.35 mM), Acrylink TFP (0.3 mM), AcX (0.35 mM)) has been reported to increase S/N in several sample types.
- Avoid contact of any dye containing reagent with the NaBH₄ solution, as many dyes will lose fluorescence.
- The NEM thiol blocking step is optional, and can be used when background staining is observed.

This product ("Product" or "Kit") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in clinical application is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to three months from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non-infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the results of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of Products.